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
**Version NO.: 1.3**

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# InstantSure Covid-19 Ag CARD

## Instruction for Use FOR IN VITRO DIAGNOSTIC USE

English

This instruction for use (IFU) must be read carefully prior to use. Instruction for use must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions for use.

### PRODUCT NAME

InstantSure Covid-19 Ag CARD

### PACKING SPECIFICATION

1 tests/ kit, 5 tests/ kit, 20 tests/ kit

### INTENDED USE

This kit is used for in vitro qualitative detection of SARS-CoV-2 antigen in human nasopharyngeal swabs, nasal swabs within 7 days after clinical symptoms.

Positive result from the test need further analyze with clinical history of patient and other diagnostic information to determine patient infection status. Positive value is only a reference guide for clinical diagnosis. The test results only reflect the current state of the sample. Negative result cannot exclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

This test is intended for use by trained clinical laboratory personnel specifically instructed and trained in the techniques of in vitro diagnostic procedures, and proper infection control procedures and individuals similarly trained in point of care settings.

### TEST PRINCIPLE

InstantSure Covid-19 Ag CARD is a colloidal gold enhanced double antibody sandwich immunoassay for the qualitative determination of SARS-CoV-2 nucleocapsid proteins antigen. Anti-SARS-CoV-2 antibody are immobilized in the test region on nitrocellulose membrane. If the specimen contains SARS-CoV-2 nucleocapsid proteins, during the assay specimen is allowed to react with the colored conjugate (SARS-CoV-2 antibody-colloidal gold conjugate); the mixture then migrates chromatographically on the membrane by the capillary action. SARS-CoV-2 positive specimen produces a distinct color band in the test region, formed by the specific antibody antigen colored conjugate complex "(Au-SARS-CoV-2-Ab) - (SARS-CoV-2-Antigen) - (SARS-CoV-2-Ab)". Absence of this colored band in the test region suggests a negative result. A colored band always appears in the control region serving as procedural control regardless of the specimen contains SARS-CoV-2 nucleocapsid proteins or not.

### REAGENTS AND MATERIALS SUPPLIED

#### 1. Main components:

Component	Specification	1 tests/kit	5 tests/kit	20 tests/kit
Test Card		1	5	20
Collection tubes include extraction solution		1	5	20
Specimen Collection Swab*		1	5	20
Package Insert		1	1	1
Work Station		The box	The box	—

\* According to customer needs

#### 2. Main ingredients of the Test card

SARS-CoV-2 antibody	Coated in the Test region on NC membrane
Goat anti Chicken IgY polyclonal antibody	Coated in the control region on NC membrane

SARS-CoV-2 antibody, Chicken IgY, Colloidal gold conjugate	Coated in the conjugate pad
Other test device supports	/

### 3. Main ingredients of the sample extraction solution

•The extraction solution is colorless or light yellow. Contains saline, detergents and preservatives.

**Note:** The components in different batches of the kit cannot be mixed.

### MATERIALS REQUIRED BUT NOT PROVIDED

• Timer •

### STORAGE AND STABILITY

- Kits should be stored at 2°C~30°C, in a cool, dark, dry place, valid for 24 months, forbidden to store under 2°C and avoid using expired products.
- The test card should be in aluminum foil bag after opening, to the specified environment (temperature 2°C~35°C, humidity 40%~60%) used within 1 hour.
- The extraction solution should be used immediately after dropping into the collection tube.
- MFD date and EXP date: marked on the label.

### SPECIMEN

- Specimens obtained early during symptom onset will contain the highest viral titers; specimens obtained after five days of symptoms are more likely to produce negative results when compared to an RT-PCR assay. Inadequate specimen collection, improper specimen handling and/or transport may yield false results; therefore, training in specimen collection is highly recommended due to the importance of specimen quality to obtain accurate test results.

Acceptable specimen type for testing is a direct swab specimen or a swab in viral transport media (VTM) without denaturing agents. Use freshly collected direct swab specimens for best test performance.

Prepare the extraction tube according to the Test Procedure and use the sterile swab provided in the kit for specimen collection.

#### 2. Nasopharyngeal Swab Specimen Collection

2.1 Remove the swab from the package.

2.2 Tilt patient's head back about 70°.

2.3 Insert the swab through the nostril parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx. (Swab should reach depth equal to distance from nostrils to outer opening of the ear.) Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions.


2.4 Slowly remove swab while rotating it. Specimens can be collected from both sides using the same swab, but it is not necessary to collect specimens from both sides if the tip of swab is saturated with fluid from the first collection. If a deviated septum or blockage creates difficulty in obtaining the specimen from one nostril, use the same swab to obtain the specimen from the other nostril.

#### 3. Nasal Swab Specimen Collection

3.1 While gently rotating the swab, insert swab about 2.5cm(1 inch) into nostril until resistance is met at turbinates.

3.2 Rotate the swab several times against nasal wall and repeat in other nostril using the same swab.

- The samples should be used as soon as possible after collected (within 1 hour).

Nasopharyngeal and Nasal specimens are stable in 30 mins when kept in the sample extraction solution provided with the kit. Samples should be tested as soon as possible after collection.


### TEST PROCEDURE

**NOTE:** Do not open the pouch until you are ready to perform a test, and the single-use test is suggested to be used under low environment humidity (RH≤70%) within 1 hour.

Allow all kit components and specimens to reach room temperature between 18°C~26°C prior to testing. Identify the test card for each specimen.

#### Direct Swab Test Procedure

- Tear off the aluminum foil film sealing the collection tube.
- Completely immerse the swab head in the extraction solution in the collection tube. Completely mix the solution by rotating the swab forcefully against the side of the tube at least 10 times (while submerged) and squeeze the tube 5 times by hand to ensure that the sample on the sampling swab is fully eluted into the extraction solution. **Leave the swab in the collection tube for one minute.**
- Squeeze the swab head along the inner wall of the collection tube to keep the liquid in the tube as much as possible. Discard the swab and cover the dropper tip to mix the liquid thoroughly.
- Remove the test card from the foil pouch and place on a clean dry surface. Dispense 80µL (3 drops) of the specimen into the circular sample well on the card.
- Interpret the test results **at 15~20 minutes**. Do not interpret the results after 20 minutes.
- Discard used test tubes and Test card in suitable biohazards waste container.


Tear off the aluminum foil film sealing the collection tube, Immerse the swab head in the collection tube, Squeeze liquid from swab, Cover the dropper tip, 80µL(3 drops) interpret the test results at 15~20 min.

#### Swab in Viral Transport Media (VTM) Test Procedure

- Insert the swab specimen into the transport tube containing a maximum of 3ml VTM without denaturing agents.

- Mix the specimen stored in VTM by vortexing.
- Transfer 300ul of the VTM solution containing specimen into the collection tube which contains extraction solution with a calibrated micropipette. Homogeneous mixture by pipetting up an down.
- Cover the collection tube with a dropper tip tightly, **and let the extracted solution stand for one minute.**

- Follow Steps 3-5 of the **Direct Swab Test Procedure** above.

### Caution:

- Do not interpret the results within 15 minutes or exceed 20 minutes.
- Use a clean pipette or tip for every sample to avoid cross-contamination.

### INTERPRETATION OF TEST RESULTS

**Positive:** Both purplish test band and purplish control band appear on the membrane.

**Negative:** Only the purplish control band appears on the membrane. The absence of a test band indicates a negative result.

**Invalid:** There should always be a purplish control band in the control region regardless of test result. If control band is not seen, it indicates that the incorrect operation process or the kit has deteriorated or damaged. In this case, read the instructions carefully again and retest with a new test card. If the problem persists, stop using product of this lot immediately and contact your local supplier.

Note: The purplish band in the test area (T) can show different color depths. However, within the specified observation time, regardless of the color of the ribbon, even a very weak ribbon should be judged as a positive result.

### LIMITATIONS

1. The result of the product should not be taken as a confirmed diagnosis, for clinical reference only. Judgment should be made along with RT-PCR results, clinical symptoms, epidemic condition and further clinical data.

1) Sensitivity of the test after the first seven days of the onset of symptoms has been demonstrated to decrease as compared to a RT-PCR SARS-CoV-2 assay.

2) Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.

3) The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection and performance may differ in asymptomatic individuals.

2. This reagent can only qualitatively detect SARS-CoV-2 antigens in human nasopharyngeal swab, nasal swab. It cannot determine the certain antigen content in viral transport media.

3.The accuracy of the test depends on the sample collection process. Improper sample collection, improper sample transportation and storage or freezing and thawing of the sample will affect the test results.

4.It is optimum when eluting swabs with the matched extraction solution. Using other diluents may result in wrong results.

5.The solution and test card must be equilibrated to room temperature (18°C~26°C) before used, otherwise the results may be incorrect.

6.Sensitivity may decrease if the sample did not test directly. Please test the sample as soon as possible.

7.Cross- reactions maybe exist due to the N protein in SARS has a high homology with the new coronavirus (SARS-CoV-2). However, the interpretation of the results is not affected during seasons without SARS infection.

8.Analysis the possibility of false negative results:

1)Inappropriate sample collection, using other non-matching solution, sample transfer time is too long (more than half an hour), the volume of solution added when eluted the swab are too much, non- standardized elution operation, low virus titer in the sample, these may all lead to false negative results.

2)Mutations in viral genes may lead to changes in antigen epitope, leading to false negative results.

9.Analysis the possibility of false positive results:

1)Inappropriate sample collection, using other non-matching solutions, non- standardized elution operation, these may all lead to false positive results.

2)Cross-contamination of samples may lead to false positive results.

3)False negative result from nucleic acid.

10.Analysis the possibility of invalid result:

1)If the sample volume is not enough, the chromatography cannot be carried out successfully.

2)The test card would invalid if the package was broken. The packaging status must be carefully checked before use.

#### PERFORMANCE CHARACTERISTIC

##### 1.Performance- bench

###### 1.1.The coincidence rate of positive controls

Tested with 5 positive controls, the results were all positive, and the coincidence rate was 5/5.

###### 1.2.The coincidence rate of negative controls

Tested with 10 negative controls, the results were all negative, and the coincidence rate was 10/10.

###### 1.3.Repeatability

Tested with repeatable control for 10 times, the results were all positive and consistent.

###### 1.4.Limit of detection

Use 3 different concentration LoD controls to test, L1 is negative, L2 ~ L3 are positive.

##### 2.Performance-Clinical

###### 2.1. For Nasopharyngeal Swab:

The clinical performance of InstantSure COVID-19 Ag CARD was established in prospective studies with nasopharyngeal swabs collected from 600 individual symptomatic patients (within 7 days of onset) and asymptomatic patients who were suspected of COVID-19.

Summary data of InstantSure COVID-19 Ag CARD as below: The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value≤33 and Ct value≤37).

COVID-19 Antigen		Clinical diagnosis PCR Ct≤33		Total
		Positive	Negative	
SAIERTECH®	Positive	264	1	265
	Negative	8	314	322
Total		272	315	587
sensitivity( Ct≤33): 97.06%(264/272), 95%CI (94.29%-98.72%) specificity (Ct≤33): 99.68%(314/315), 95%CI(98.24%-99.99%)				

COVID-19 Antigen		Clinical diagnosis PCR Ct≤37		Total
		Positive	Negative	
SAIERTECH®	Positive	267	1	268
	Negative	18	314	332
Total		285	315	600
sensitivity(Ct≤37): 93.68%(267/285), 95%CI (90.20%-96.21%) specificity(Ct≤37): 99.68% (314/315), 95%CI(98.24%-99.99%)				

###### 2.2 For Nasal Swab:

The clinical performance of InstantSure COVID-19 Ag CARD was established in prospective studies with nasal swabs collected from 600 individual symptomatic patients (within 7 days of onset) and asymptomatic patients who were suspected of COVID-19.

Summary data of InstantSure COVID-19 Ag CARD as below: The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value≤33 and Ct value≤37).

COVID-19 Antigen		Clinical diagnosis PCR Ct≤33		Total
		Positive	Negative	
SAIERTECH®	Positive	262	2	264
	Negative	10	313	323
Total		272	315	587
sensitivity( Ct≤33): 96.32%(262/272), 95%CI (93.34%-98.22%) specificity(Ct≤33): 99.37%(313/315), 95%CI(97.73%-99.92%)				

COVID-19 Antigen		Clinical diagnosis PCR Ct≤37		Total
		Positive	Negative	
SAIERTECH®	Positive	265	2	267
	Negative	20	313	333
Total		285	315	600
sensitivity( Ct≤37): 92.98%(265/285), 95%CI (89.37%-95.66%) specificity(Ct≤37): 99.37%(313/315), 95%CI(97.73%-99.92%)				

###### 2.3 For all samples including nasopharyngeal and nasal swabs:

Summary data of InstantSure COVID-19 Ag CARD as below:

COVID-19 Antigen		Clinical diagnosis PCR Ct≤33		Total
		Positive	Negative	
SAIERTECH®	Positive	526	3	529
	Negative	18	627	645
Total		544	630	1174
sensitivity( Ct≤33): 96.69%(526/544), 95%CI (94.82%-98.03%) specificity(Ct≤33): 99.52%(627/630), 95%CI(98.61%-99.90%)				

COVID-19 Antigen		Clinical diagnosis PCR Ct≤37		Total
		Positive	Negative	
SAIERTECH®	Positive	532	3	535
	Negative	38	627	665
Total		570	630	1200
sensitivity( Ct≤37): 93.33%(532/570), 95%CI (90.96%-95.24%) specificity(Ct≤37): 99.52%(627/630), 95%CI(98.61%-99.90%)				

##### 3. Cross-reactivity

Cross-reactivity of the cassette was evaluated. The results showed no cross reactivity with the following microorganism.

No.	Virus name	Test conc.
1	Human coronavirus 229E	1.6×10 <sup>5</sup> TCID <sub>50</sub> /mL

2	Human coronavirus OC43	1.6×10 <sup>4</sup> TCID <sub>50</sub> /mL
3	Human coronavirus NL63	1.6×10 <sup>3</sup> TCID <sub>50</sub> /mL
4	MERS-coronavirus (Irradiated)	8.9×10 <sup>4</sup> TCID <sub>50</sub> /mL
5	Adenovirus	1.6×10 <sup>5</sup> TCID <sub>50</sub> /mL
6	Human metapneumovirus (hMPV)	1.6×10 <sup>5</sup> TCID <sub>50</sub> /mL
7	Human parainfluenza virus 1	8.9×10 <sup>5</sup> TCID <sub>50</sub> /mL
8	Human parainfluenza virus 2	1.0×10 <sup>5</sup> TCID <sub>50</sub> /mL
9	Human parainfluenza virus 3	1.6×10 <sup>5</sup> TCID <sub>50</sub> /mL
10	Human parainfluenza virus 4 a	1.6×10 <sup>3</sup> TCID <sub>50</sub> /mL
11	Human parainfluenza virus 4 b	5.0×10 <sup>5</sup> TCID <sub>50</sub> /mL
12	Influenza A	5.2×10 <sup>5</sup> TCID <sub>50</sub> /mL
13	Influenza B	1.8×10 <sup>5</sup> TCID <sub>50</sub> /mL
14	Enterovirus	1.6×10 <sup>5</sup> TCID <sub>50</sub> /mL
15	Respiratory Syncytial Virus A	2.8×10 <sup>5</sup> TCID <sub>50</sub> /mL
16	Rhinovirus 16	5×10 <sup>3</sup> TCID <sub>50</sub> /mL
17	Haemophilus influenzae	800 cfu/vial
18	Streptococcus pneumoniae	4.25×10 <sup>5</sup> CFU/mL
19	Streptococcus pyogenes	800 cfu/vial
20	Bordetella pertussis	4.8×10 <sup>6</sup> CFU/mL
21	Mycoplasma pneumoniae	3.0×10 <sup>5</sup> CFU/mL
22	Chlamydia pneumoniae	9.1×10 <sup>6</sup> IFU/mL
23	Legionella pneumophila	3.9×10 <sup>5</sup> CFU/mL
24	Staphylococcus aureus	800 cfu/vial
25	Staphylococcus epidermidis	800 cfu/vial
26	Candida albicans	N/A

##### 4.Interference Substances

The test results do not be interfered with the substance at the following concentration:

No.	Substances	Conc.
1	Whole Blood	4%
2	Chloraseptic(Menthol/Benzocaine)	1.5mg/mL
3	Naso GEL(NeiMed)	5% v/v
4	CVS Nasal Drops (Phenylephrine)	15% v/v
5	Afrin (Oxymetazoline)	15% v/v
6	CVS Nasal Spray (Cromolyn)	15% v/v
7	Zicam	5% v/v
8	Homeopathic (Alkalol)	1:10 dilution
9	Sore Throat Phenol Spray	15% v/v
10	Tobramycin	4 µg/mL
11	Mupirocin	10 mg/mL
12	Fluticasone Propionate	5% v/v
13	Tamiflu (Oseltamivir Phosphate)	5 mg/mL
14	Mucin	0.5%

##### PRECAUTIONS

- The reagent is a disposable diagnostic reagent for in vitro.
- The operation should be carried out strictly according to the instructions. Do not use expired or damaged products.
- Reagents should be used as soon as possible (within 1 hour) after removal from aluminum foil bags, so as to avoid exposure to air for too long and affecting test results due to dampness.
- Do not use samples that have been placed for too long or contaminated.
- Please operate in accordance with the laboratory testing procedures for infectious diseases. Waste after use should be treated in accordance with infectious substances and should not be discarded at will.
- Incorrect operation may affect the accuracy of the results, such as insufficient sample mixing, insufficient amount, inaccurate detection time, etc.
- Components in different batch should not be mixed.

8. There should be appropriate biosafety assurance procedures for those substances containing and suspected sources of infection. The following are relevant considerations:

- Handle samples and reagents with gloves;
- Do not suck samples with your mouth;
- Do not smoke, eat, drink, cosmetic or handle contact lenses while handling these items;
- Disinfect the spilled sample or reagent with disinfectant;
- Disinfect and treat all samples, reagents and potential pollutants in accordance with relevant local regulations;
- Each component of the reagent remains stable until the expiry date under proper handling and storage conditions. Do not use the expired reagent kit.

##### SYMBOLS


Symbol	Used for	Symbol	Used for
	Use-by date		Consult instructions for use
	Batch code		In vitro diagnostic medical device
	Temperature limit		Manufacturer
	Tests per kit		Authorized representative in the European Community
	Please don't reuse it		Don't use the product when the package is damaged

##### 【BASIC INFORMATION】

**Suzhou Soochow University Saier Immuno Biotech Co., Ltd.**

No.502 Zhongnan Road, Suzhou Industrial Park,215123, Suzhou, Jiangsu, China.

**Lotus NL B.V.**

Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

Email: peter@lotusnl.com

##### 【INSTRUCTION APPROVAL】

Version NO.: 1.0

Effective Date: Aug 11, 2021